- ♦ **The Wallaby® 3 Phototherapy System
Professional Manual**
- ♦ **Système de photothérapie Wallaby® 3
Manuel du revendeur**
- ♦ **Wallaby® 3 Phototherapiesystem
Gebrauchsanweisung**
- ♦ **El sistema de fototerapia Wallaby® 3
Manual profesional**

3606

CE
0459

Warning

This guide is not meant as a substitute for the services of a physician or other healthcare professional. The Wallaby® 3 Phototherapy System was developed to provide effective, state of the art phototherapy as a safe, comfortable treatment method of jaundice in newborn babies.

Warranty

Respironics warrants your Wallaby® Phototherapy System against defects in material and workmanship of the Illuminator unit only for a period of one (1) year from the date of purchase, one hundred eighty (180) days on the fiberoptic panels, EG-2000 and EG-2000N. This warranty does not cover the lamp, or any damage to the illuminating unit or the fiberoptic panel caused by accident, misuse, tampering, or negligence such as failure to follow the instructions provided in this guide. In the event your phototherapy illumination unit fails to give satisfactory performance within the warranty period and conditions, Respironics Georgia will repair or replace your illuminating unit at no charge for parts or labor. The foregoing warranties are in lieu of all other warranties expressed or implied, including without limitation any implied warranty of merchantability or fitness for a particular purpose.

Respironics Georgia, Inc. 1255 Kennestone Circle Marietta, GA 30066 Phone: 800-345-6443	Respironics Deutschland Gewerbestrasse 17 82211 Hersching Germany +49 8152 93060	Respironics Asia Pacific 3/F Microtron Building 38 Hung To Road Kwun Tong, Kowloon, Hong Kong
--	---	--

Product Responsibility

This Product will perform in conformity with the description thereof contained in this manual and accompanying labels and/or inserts, when assembled, operated, maintained and repaired in accordance with the instructions provided. This Product must be checked periodically. A defective Product should not be used. Parts that are broken, missing, plainly worn, distorted or contaminated should be replaced immediately. Failure to use Respironics replacement parts, including lamps, renders any written or implied warranty on the unit null and void. Should repair or replacement become necessary, Respironics recommends that a telephone or written request for service advice be made to the nearest Respironics Service Center. This Product or any of its parts should not be repaired other than in accordance with written instructions provided by Respironics and by Respironics trained personnel. The Product must not be altered without the prior written approval of Respironics Quality Assurance Department. The user of this Product shall have the sole responsibility for any malfunction that results from improper use, faulty maintenance, improper repair, damage, or alteration by anyone other than Respironics.

Wallaby® is a registered trademark of Respironics Inc.

Table of Contents

About this Manual	1
Before Using the Wallaby Phototherapy System	1
Overview	3
What is Jaundice?	3
What is Bilirubin?	3
What is Phototherapy?	3
What is the Wallaby 3 Phototherapy System?	3
Warnings and Cautions	5
Warnings	5
Cautions	6
Symbols	6
Electromagnetic Compatibility	6
Classification	6
Setup	7
Checking out the System	7
Setting Up the Wallaby System	7
Patient Setup - Wrap-Around Panel	9
Patient Setup - Flat Neonatal Panel	9
Maintenance	11
Maintenance Schedule	11
Cleaning the Illuminator Unit and Fiber Optic Panel	11
Cleaning the Optical Filter	12
Replacing the Lamp	12
Preserving Lamp Life	13
Servicing	15
Equipment Servicing	15
Troubleshooting	17
Appendix A: Joey™ Dosimeter	19
What is the Joey Dosimeter?	19
Checking Light Intensity	19
Obtaining an Intensity Reading	20
For Flat Light Panel EG-2000N	20
For Wrap Around Panel EG-2000	20
Maintenance	21
Digital Display	21
Battery	21
Recalibration	21
Wallaby 3 Specifications	23
Ordering Information	25
Index	27

Table of Contents

About this Manual

This manual is written for the medical equipment dealer and describes how to operate, maintain, and service the Wallaby 3 Phototherapy System.

Before Using the Wallaby Phototherapy System

- Read through sections 1 through 3 of this manual. Pay special attention to the **Warnings and Cautions** that appear.
- Read the **Product Responsibility** statement located on the inside front cover. It describes what is expected of you to maintain a safe and effective product.
- Read the **Warranty**. **It describes Respiromics' responsibility in case of a functional defect.**

If you have any questions about the product, contact your Respiromics Service Representative.

About this Manual

Overview

This section explains how the Wallaby Phototherapy System is used to treat infant jaundice.

What is Jaundice?

Jaundice is a yellowish skin discoloration caused by an excess of bilirubin that has accumulated in the fatty tissues of the skin.

What is Bilirubin?

Bilirubin is a naturally occurring molecule, released into the bloodstream when red blood cells break down. The breaking down of red blood cells occurs as a normal part of the constant renewing of our bodies. Usually, the level of bilirubin is of no concern since the liver breaks down the bilirubin so that it is excreted through urine and bile.

Infants are born with extra red blood cells that provide the oxygen needed while the baby grows in the mother's uterus. Once the baby is born, the excess red blood cells break down and bilirubin begins to build up. Some newborns have a premature liver that cannot effectively process the extra bilirubin, jaundice results. It is normal in all newborns for bilirubin to rise during the first 3 to 4 days of life, peak, and slowly decline. In a few babies, the bilirubin rises beyond what is considered an acceptable level and therapy is required.

Left untreated, high bilirubin levels can be toxic and can possibly cause brain damage. Phototherapy should be initiated when the infant's serum bilirubin concentration exceeds the acceptable limits determined by the baby's physician. After treatment has begun, the baby's bilirubin level should be measured to make sure that the therapy is effective. The bilirubin should drop to an acceptable level in a few days.

What is Phototherapy?

Phototherapy simply means treatment with light. In phototherapy, the special light causes the bilirubin to break down into a form that can be excreted by the baby.

Phototherapy previously required the baby be readmitted to the hospital. Today, with the development of portable lights and systems such as the Wallaby Phototherapy System, babies can be kept at home while treatment takes place.

What is the Wallaby 3 Phototherapy System?

Any phototherapy system uses visible light to convert bilirubin to waste products that are mostly excreted into bile, thus reducing the bilirubin level in the baby's blood.

The Wallaby is a fiberoptic phototherapy system that consists of an Illuminator unit – the source of the therapeutic light – and a fiberoptic panel that delivers the therapeutic light to the baby.

The Illuminator unit sends light out the fiberoptic cable to the entire area of the panel. The panel is inserted into a protective cover composed of hypoallergenic material or soft cotton. This wrap is soft and comfortable and allows the therapeutic light to be emitted towards the baby.

Overview

With the design of the Wallaby System, the baby can be held and fed and enjoy the healing comfort of parents and home while treatment is administered. Additionally, when the Wallaby system is properly used the baby's eyes do not need to be protected as with conventional phototherapy.

Intended Use:

The Wallaby 3 is intended to treat hyperbilirubinemia through phototherapy in the home or hospital.

Warnings and Cautions

Warnings

A *warning* indicates a potentially harmful situation.

Federal law (USA) restricts this device to sale by, or on the order of, a licensed physician. This instrument should be used only under the supervision of a physician.

1. Use the Wallaby only for its intended use as described in this manual.
2. Eye patches are only necessary if the panel will be illuminated while not attached to the patient. Normal use of the Wallaby Illuminator with the panel in a disposable or reusable cover does not require the use of eye patches.
3. Do not leave the Wallaby unit on when the fiberoptic panel is not around the baby.
4. Always turn off and unplug the Wallaby unit when cleaning, changing the lamp, or servicing the unit.
5. Do not use the Wallaby unit in the presence of flammable substances such as anesthetics or cleaning agents.
6. Do not use the Wallaby unit where oxygen is administered.
7. Do not use while bathing the baby.
8. Do not place or store the Wallaby unit where it can fall or be pulled into a tub or sink.
9. If the Wallaby unit falls into water or if fluid is spilled on the device do not reach for it without first unplugging the cord. Discontinue use of the device and contact an authorized Respiration service center.
10. Always connect the Wallaby unit to a properly grounded outlet.
11. Never operate the Wallaby if it has a damaged cord or plug. Do not insert anything into the end of the plug.
12. Do not use an extension cord.
13. Do not walk around with the baby wrapped in the fiberoptic panel.
14. Carefully place the panel cord to avoid entanglement.
15. Always place the Illuminator on a hard flat surface.
16. Do not place the Illuminator unit in an isolette or infant warmer. Do not place a temperature sensor for these devices under the panel.
17. Do not operate the Wallaby unit on its side.
18. Never block the air vents of the unit or place it on a soft surface such as a bed, crib, carpeted floor, or couch where the air vents may be blocked.
19. Never put fingers or objects into the opening on the Wallaby unit.
20. The inside of the Wallaby unit becomes very hot during use. Wait 10 minutes before changing light bulbs.
21. The Wallaby device utilizes high voltage, refer servicing to qualified personnel.
22. The reusable panel covers are for use in home care phototherapy only.

Cautions

A *caution* indicates a condition that may lead to equipment damage, malfunction, or inaccurate operation.

1. Keep the unit away from any heated surface.
2. Do not scratch, damage, or soil the ferrule end of the panel. Also, do not place sharp or heavy items on the panel, this can damage the panel and affect its light output.
3. Do not dry the fiberoptic panel with artificial heat.
4. Allow the unit to cool 10 minutes before moving.

Symbols

	Attention: read accompanying documents
	Type BF applied part
	Caution, Hot Surface
	Caution, High Voltage
	Alternating Current
REF	Part Number
SN	Serial Number
I	On (power:connection to the mains)
O	Off (power: connection to the mains)

Electromagnetic Compatibility

The Wallaby 3 Phototherapy System complies with Electromagnetic Compatibility Standards for Medical Devices. If you suspect that the Wallaby 3 is interfering with the operation of other appliances or if they are interfering with the Wallaby 3:

- Move the Wallaby 3 away from the appliance
- Connect the Wallaby 3 and the appliance to different AC outlets so that they are operating off of different circuit breakers

Classification

The Wallaby 3 Phototherapy System is classified as:

- IEC Class I Equipment
- Type BF Applied Part
- Ordinary protection against the ingress of liquids
- Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide
- Continuous operation

Setup

Checking out the System

Prior to each patient setup, the light intensity of the unit should be checked. The Joey Dosimeter should be used to measure the light intensity. For more information on using the Joey Dosimeter refer to Appendix A.

Setting Up the Wallaby System

1. Place the Illuminator unit on its bottom on a hard, flat surface, away from any heat source. Make sure it is no more than four feet from where the baby will be positioned.



IMPORTANT! Do not block any of the air vents on the sides of the Illuminator unit.


IMPORTANT! Do not place the Illuminator unit on its side.

2. Insert the metal end of the cable, called the ferrule, with the metal post facing up, into the opening. Twist the panel $\frac{1}{4}$ turn to lock into place. The light will not come on if the panel is not inserted in the Illuminator.
3. Plug the Illuminator unit into an electrical outlet.
4. Place the panel on the baby. (See Patient Setup-Wrap Around Panel or Neonatal Panel.)
5. Be sure the Light Level on the Illuminator unit is at the setting recommended by the physician. Then, turn the Illuminator unit ON.

- 6.** After the power is turned on, the green power light will illuminate. The System Alert light will indicate the bulb status.

- If one bulb is burned out, the yellow light will flash. (See Troubleshooting.)
- If both bulbs have burned out or if the unit overheats, the yellow light is solid. (See Troubleshooting.)

IMPORTANT! If the power cord needs repair or replacement, do not connect the unit.

IMPORTANT! Connect the power cord to a properly grounded outlet. Make sure that the plug contains the ground pin.

Patient Setup - Wrap-Around Panel


This section explains how to prepare a baby for a phototherapy treatment using the wrap-around fiberoptic light panel.

This panel provides full coverage. It is most commonly used on patients at home.

1. Insert the panel into a disposable or single-patient use cover with the light facing the sheer side of the cover.
2. Put a T-shirt on the baby and roll it up from the bottom until it fits comfortably under the baby's arms.
3. Place the covered panel under the baby's torso, positioning it so it is under the baby's armpits.
4. Wrap the panel around the baby as shown in the illustration. Use the tape or Velcro tabs to secure the panel around the baby.
5. If the cover becomes soiled, discard it (disposable) or replace it with a clean one and wash the soiled cover (cloth).

- For a larger or more active baby, you may want to tape the panel to the baby's diaper.
- To be sure the panel is not wrapped too tightly, insert your finger between the panel and the baby's body. Your finger should fit easily.
- You may wrap your baby in a blanket or put the baby in a sleeper.

Patient Setup - Flat Neonatal Panel

This section explains how to prepare a baby for a phototherapy treatment using the flat neonatal fiberoptic light panel.

This panel is ideal for preemie or underweight infants; it may also be used on full term infants.

1. The protective cover for the neonatal panel is a T-vest. Insert the panel in the vertical section of the T with the light facing the sheer side of the cover.
2. Lay the covered panel on a flat surface. Be sure the light emitting side is facing up.
3. Position the baby's chest or back directly on the panel. The cable connected to the panel should be between the baby's legs.
4. Secure the T-vest to the baby by first wrapping the side without the tape around the baby's midsection. Then, wrap the side with the tape tab over the infant, peel off the protective cover on the tab, and secure it. Be sure the vest is snug.

5. If the T-vest becomes soiled, discard it and replace it with a new one.

- For a larger or more active baby, you may want to tape the panel to the baby's diaper.
- Make sure the panel is not wrapped too tightly, insert your finger between the panel and the baby's body. Your finger should fit easily.
- You may wrap your baby in a blanket or put the baby in a sleeper.

Maintenance

This section explains how to clean, check, and maintain the Wallaby system.

Maintenance Schedule

Follow the schedule and instructions below to keep the Wallaby system properly maintained.

After each patient	As needed
<ol style="list-style-type: none"> 1. Check light intensity (See Appendix A). 2. Replace the disposable light panel cover or vest. 3. Inspect the power cord and plug for damage. 4. Verify light doesn't illuminate when panel is not inserted. 	<ol style="list-style-type: none"> 1. Clean the Illuminator unit and fiber optic panel. 2. Replace the lamp. 3. Clean the optical filter.

IMPORTANT! An authorized service representative should replace all parts, except the lamp.

Cleaning the Illuminator Unit and Fiber Optic Panel

Follow the instructions in this section any time the Illuminator unit or fiber optic panel is dirty.

IMPORTANT! Be sure the Illuminator unit is turned off and is unplugged before cleaning. Do not immerse any part of the equipment in any liquid.

- IMPORTANT!** Keep panel away from sharp objects that could scratch or puncture the cover.
- Use soapy water, a 10% bleach solution or full strength ammonia.
 - Use a soft sponge or cloth to apply the cleaner.
 - Apply the product to the sponge or cloth and wipe down the fiberoptic panel and Illuminator.
 - Allow the equipment to air dry. DO NOT DRY WITH ANY MEANS OF ARTIFICIAL HEAT.

DO NOT USE

- Phenolic compound based germicide cleaner/disinfectant
- Gluteraldehyde disinfectant/sterilants
- Regular commercial cleaners or laundry detergents
- Iodine solutions, strong acids or strong alkali solutions

These solutions could leave residue on the surfaces, and /or be abrasive or harmful to the infant.

Cleaning the Optical Filter

The optical filter is located inside the Illuminator unit. You should clean the mirror whenever your irradiance levels are below the acceptable range.

1. Using a fresh cotton swab, wet with alcohol, wipe the mirror in one direction.
2. After one pass across the face of the mirror, throw the swab away. DO NOT scrub from side to side.
3. Use a fresh swab and repeat Steps 1 and 2 on the other side of the mirror.

IMPORTANT! Be sure the Illuminator unit is turned off and is unplugged before cleaning. Do not immerse any part of the equipment in any liquid.

IMPORTANT! Do not touch the optical filter with your hands. Oil from skin can degrade the optical filter.

Replacing the Lamp

Part Number 100-3600-20

The Wallaby is equipped with a System Alert light that blinks whenever a lamp fails to operate. Two lamps are built in. When one fails, the second lamp automatically switches into place.

IMPORTANT! Replace the lamps ONLY with the lamp specified. Use of any other lamp than the one specified in this manual will affect the performance of the system and could cause damage to the Illuminator or fiber optic panel. Such use could also decrease the effectiveness and safety of therapy.


IMPORTANT! Failure to use Resironics, Inc. manufactured lamps will render the warranty on the Wallaby Phototherapy System null and void.

To replace a burned out bulb

1. Turn off the Illuminator unit and disconnect the power cord from the power source. Wait 5 minutes for lamp to cool, if appropriate.
2. Remove the 2 screws and remove lamp access panel. A torx head screwdriver size 15 is required.
3. Remove the burned out bulb by pulling the light straight out. This is the bulb either on the far right or far left of the system. The good bulb is in the center.
4. Replace with a new light bulb (P/N 100-3600-20).
5. Replace the access panel and screws.

There is no warm up time required with a new bulb as is the case with conventional phototherapy. Therapy is immediately effective.

Preserving Lamp Life

Read this section for guidelines in preserving lamp life. Proper light output directly impacts the effectiveness of treatment.

Factors, which affect light output and lamp life:

- Light intensity setting
- Vibration and mechanical shock
- Soil on lamp and/or heat reflector

Light Intensity Selector Setting

Lamp life will vary with the intensity at which the Illuminator unit is operated. Approximate lamp life:

- Level 1 = 600 hours/bulb
- Level 2 = 285 hours/bulb

Vibration and Mechanical Shock

To maximize the life of the bulb, avoid placing the Illuminator unit on a surface that is not stable.

Soil on Lamp

Using a fresh cotton swab, wet with alcohol, swab the lamp. Allow the lamp to air dry before turning the Illuminator on. Do not touch or leave fingerprints on lamp bulb. Oil from the skin can reduce lamp life.

IMPORTANT! Be sure the Illuminator unit is turned off and is unplugged before cleaning. Do not immerse any part of the equipment in any liquid.

Maintenance

Servicing

Equipment Servicing

Repair and service of equipment under warranty should be performed at the Respiromics, Inc Service Center. Service performed or attempted by unauthorized personnel may void the warranty. To promote full reliability and ensure performance to factory specifications, have all repairs and service performed by authorized Respiromics service representatives. For Technical Support, call 800-421-8754.

Servicing

Troubleshooting

Following is a list of problems that may occur while using the Wallaby 3 System. For additional information, please contact your Resironics Technical Support representative.

Problem	Reason/Action
Power Light is not green	<ul style="list-style-type: none"> <input type="checkbox"/> Check to make sure power cord is properly attached and plugged into an active electrical outlet. <input type="checkbox"/> Ensure the unit is turned on.
System Alert LED is blinking yellow	<ul style="list-style-type: none"> <input type="checkbox"/> One light bulb has burned out and the system has automatically switched to the second light bulb. Continue therapy until completion and replace lamp before next patient.
System Alert LED is solid yellow	<ul style="list-style-type: none"> <input type="checkbox"/> Unit may have overheated or both light bulbs have burned out. Allow unit to cool, if condition persists replace light bulbs.
Light is not being emitted from fiber optic panel	<ul style="list-style-type: none"> <input type="checkbox"/> Check to make sure panel is securely locked into Illuminator box. <input type="checkbox"/> Check System Alert LED. If solid yellow, both lights are burned out or unit has overheated. Allow unit to cool/replace light bulbs.
Light is emitted from Illuminator unit when the fiberoptic panel is not inserted	<ul style="list-style-type: none"> <input type="checkbox"/> Replace the panel input disconnect inside the unit.

Troubleshooting

Appendix A: Joey™ Dosimeter

What is the Joey Dosimeter?

The Joey Dosimeter is a wide-band *spectroradiometer* used to measure the light wavelengths emitted from the fiber optic light panel in a phototherapy system. Checking light output helps ensure proper results from therapy.

The meter measures light intensity in units of microwatts per square centimeter per nanometer ($\mu\text{W}/\text{cm}^2/\text{nm}$). The Joey has one scale ranging from 00.0 to 99.9 $\mu\text{W}/\text{cm}^2/\text{nm}$ (425-550 nanometers). The light intensity ($\mu\text{W}/\text{cm}^2/\text{nm}$) is read directly from the digital display.

You should measure light output whenever there is a change that affects the phototherapy light output and prior to each new patient use. The light output should be checked on both Light Level 1 and Light Level 2. See Checking Light Intensity to verify minimum acceptable light level.

Joey™ is a trademark of Respiromics, Inc.

Checking Light Intensity

You should measure the light output of the Wallaby System:

- Prior to each new patient use
- After anything occurs that could impact light output, such as
 - The cleanliness of the fiber optic light panel is questionable
 - The fiberoptic light panel has been scratched or otherwise damaged
 - There has been a change in the line voltage

IMPORTANT! Variations within the lamps can affect the lamp life and light output. The vibrations and mechanical shocks can also affect the lamp life. Do not use non-recommended lamps at any time. These lamps may not only affect the light output, but may also affect the entire phototherapy system.

IMPORTANT! To ensure accurate light intensity measurements the Joey Dosimeter must be calibrated yearly.

When checking the light output, be sure to:

- Check light output with the Illuminator unit set at both light intensity levels:**
Level 1 and Level 2. The minimum readings should be:

	Level 1	Level 2
Standard Wrap Around Panel EG-2000	10 μ W/cm ² /nm*	15 μ W/cm ² /nm*
Neonatal Flat Panel EG-2000N	25 μ W/cm ² /nm*	35 μ W/cm ² /nm*

*Minimum average irradiance level of the illuminated area.

Obtaining an Intensity Reading

Follow the steps listed below to obtain a light reading using the Joey Dosimeter.

For Flat Light Panel EG-2000N

1. Be sure the light panel and the cable are properly attached to the Illuminator.
2. Take the protective cover off of the fiberoptic panel.
3. Plug in the Illuminator unit.
4. Turn the Illuminator unit ON.
5. Set Light Level to 1.
6. Place sensor of Joey face down on one end of the light emitting side of panel.
7. Depress the PUSH TO READ button on the Joey.
8. Note the reading on the Joey display.
9. Repeat steps 5-7 as 3 times across the panel.
10. Average the 3 readings together and ensure the value is at least 25 μ W/cm²/nm.
11. Repeat the test on Light Level 2 and ensure the value is at least 35 μ W/cm²/nm.

For Wrap Around Panel EG-2000

1. Be sure the light panel and the cable are properly attached to the Illuminator unit.
2. Take the protective cover off of the fiberoptic panel.
3. Plug in the Illuminator unit.
4. Turn the Illuminator unit ON.
5. Set the Light Level to 1.
6. Place sensor of Joey face down on one end of light emitting side of panel.
7. Depress the PUSH TO READ button on the Joey.
8. Note the reading on the Joey display.
9. Repeat steps 5-7 3 times across the panel.
10. Average the 3 readings together and ensure the value is at least 10 μ W/cm²/nm.
11. Repeat the test on Light Level 2 and ensure the value is at least 15 μ W/cm²/nm.

Maintenance

Follow the instructions in this section to ensure proper maintenance of the Joey Dosimeter.

Digital Display

- Avoid getting fingerprints, dirt, grease or dust on the photosensor.
- A dirty sensor can result in false low intensity readings.
- Clean the sensor with a water-dampened lens paper or soft cloth. Then, wipe dry.

WARNINGS and CAUTIONS

- DO NOT immerse the photosensor in water or any other liquid.
- DO NOT apply water or commercial cleaning agents directly to the sensor and/or digital display.
- Use a soft damp cloth or lens paper to wipe the display window while exerting very little pressure.
- DO NOT rub your fingers across the display window. This creates static electricity and can cause false readings.

Battery

The Joey Dosimeter operates with a single 9-volt, non-rechargeable battery. A LO BAT signal will display if the battery needs to be changed.

Follow the steps shown below to replace the battery:

1. Remove the battery compartment cover from the lower backside of the Joey Dosimeter.
2. Lift out the 9-volt battery.
3. Replace it with a battery of the same type.
4. Replace the battery compartment cover.

Recalibration

The Joey Dosimeter is calibrated against lamp standards traceable to the National Institute of Standards Technology. It is calibrated to be +/- 2% or +/- one digit (the least significant digit).

Proper calibration can only be accomplished with a National Institute or Standards Technology traceable lamp standard. It is recommended that the instrument be returned to Respiration for recalibration immediately if the meter has been physically abused.

Recommended frequency for recalibration is annually.

To have the Joey Dosimeter serviced or calibrated:

- Call Customer Service at 800-345-6443 for a return authorization number.
- Securely pack the dosimeter and put the return authorization number on the outside of the package.

Wallaby 3 Specifications

Illuminator:

Model: 3600

Size: 10.25"W x 9.0"D x 4.0"H (26.67cm W x 22.86cm D x 10.16cm H)

Weight: 7.15 lb. (3.24 kg)

Fiberoptic Panel:

Model: EG-2000 (Standard Panel)

Overall Pad Size-Standard: 4" x 15" (excluding boot) (10.16cm x 38.10cm)

Illuminated Area-Standard: 3" x 14" (7.62cm x 35.56cm)

Model: EG-2000N (Neonatal Panel)

Overall Pad Size-Neonatal: 5" x 7" (excluding boot) (12.70cm x 17.78cm)

Illuminated Area-Neonatal: 4" x 6" (10.16cm x 15.24cm)

Irradiance Level:

Standard Panel-Light Level 1: 10 μ W/cm²/nm*

Standard Panel-Light Level 2: 15 μ W/cm²/nm*

Neonatal Panel-Light Level 1: 25 μ W/cm²/nm*

Neonatal Panel-Light Level 2: 35 μ W/cm²/nm*

*Minimum acceptable irradiance level.

Storage Temperature:

-20 to +60 degrees C, RH 95%

Operating Temperature:

15 to 35 degrees C, RH 95%

Electrical Requirements:

Nominal Input Voltage: 100-230 VAC, 50/60 Hz

Nominal Power Consumption: 200 watts

Ordering Information

Description	Part Number
Wallaby 3 Phototherapy System	3600
Panel, Wrap Around	EG-2000
Panel, Neonatal	EG-2000N
Joey Dosimeter	JD-101
Disposable Covers, Wrap-Around for EG-2000	MD-510
Disposable Covers, Vest for EG-2000N	MD-610
Parents' Guide	PGM302
Replacement Phototherapy Light Bulb	100-3600-20
Replacement Power Cord	161-5000-00

Ordering Information

Index

B

Bilirubin, 3

C

Caution, 6

E

Electromagnetic Compatibility Standards, 6
Equipment, cleaning, 13
Equipment, repair, 17
Equipment, servicing, 17

F

Fiberoptic panel, 3

I

Illuminator unit, 3

J

Jaundice, ii, 3
Joey Dosimeter, 9, 21, 22, 23
Joey Dosimeter, maintenance, 23
Joey Dosimeter, recalibration, 23

L

Lamp life, 15
Lamp replacement, 14
Light intensity, 13
Light intensity reading, wrap-around panel, 22
Light intensity, measuring, 9
Light output, measuring, 21
Light wavelengths, measuring, 21

M

Maintenance schedule, 13

P

Patient setup, 9
Patient setup - flat neonatal panel, 11
Patient setup - wrap around panel, 11
Phototherapy, definition, 3
Product responsibility, ii

R

Recalibration, light meter, 23

S

Symbols, 6

T

Troubleshooting, 19

W

Wallaby 3, equipment classification, 6
Wallaby 3, intended use, 4
Wallaby Phototherapy System, 3
Wallaby Specifications
 Electrical Requirements, 25
 Fiberoptic Panel, 25
 Illuminator, 25
 Irradiance Level, 25
 Operating Temperature, 25
 Storage Temperature, 25
Warning, 5
Warranty, ii, iii, 1

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